

Distributed By:



Tel: 1.800.258.3476 | Fax: 1.866.216.4161

AXIS® TUTOPLAST® PROCESSED DERMIS



Read this entire package insert carefully prior to use.



Single patient use only, on a single occasion.



Restricted to sale by or on the order of a physician.

DESCRIPTION

Axis Tutoplast Processed Dermis is dehydrated, *Tutoplast* processed dermis from donated human tissue. The implant is preserved by the Tutoplast® process which retains the three-dimensional collagen structure responsible for the omni-directional, mechanical properties of the original dermal tissue.

Axis Tutoplast Processed Dermis is restricted to homologous use for the repair, replacement, reconstruction or augmentation of soft tissue by a qualified healthcare professional (i.e., physician). This includes supplemental support and reinforcement of soft tissue and grafting for horizontal and vertical soft tissue augmentation of thickness and length. The implant is provided sterile and requires rehydration prior to use.

DONOR SCREENING AND TESTING



This symbol on the outer label indicates the unique number assigned to the tissue donor.

The donated human tissue utilized for this implant was recovered from a donor screened for risk factors associated with infectious diseases and medical conditions that rule out donation. The donor's blood was tested for relevant communicable diseases in a laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 or equivalent, and registered with the US Food and Drug Administration for donor testing. The following test criteria were met for this donor:

REQUIRED INFECTIOUS DISEASE TESTING	
BLOOD TEST	ACCEPTABLE RESULT
HIV-1/ HIV-2 Antibody	Negative/ Non-Reactive
Hepatitis C Virus Antibody	Negative/ Non-Reactive
Hepatitis B Surface Antigen	Negative/ Non-Reactive
Hepatitis B Core Antibody (Total)	Negative/ Non-Reactive
Syphilis	Negative/ Non-Reactive
Human T-Cell Lymphotropic Virus I / II Antibody	Negative/ Non-Reactive
HIV-1 NAT-TMA	Negative/ Non-Reactive
HCV NAT-TMA	Negative/ Non-Reactive

A licensed physician for RTI Biologics, Inc. determined that the donor met eligibility requirements. The physician utilized available relevant information which may have included, but was not limited to: family/next-of-kin interview, medical/hospital records, donor physical assessment, infectious disease test results, radiology/ pathology reports, death certificate and autopsy report (if performed).

WARRANTY STATEMENT

*This biologic graft, processed and packaged for surgical implantation, is unique and does not constitute a product under liability laws. No implied warranties of merchantability or fitness for a particular purpose are applicable. No implied warranties exist as to defects in biologics which cannot be detected, removed, or prevented by reasonable use of available scientific procedures or techniques. Furthermore, **ALL WARRANTIES ARE DISCLAIMED, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE INCLUDING ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.***

PROCESSING



This symbol on the outer label indicates a unique serial number assigned to the implant, used for traceability.

The implant was processed in a controlled environment from a single donor. Microbial testing was performed where appropriate and results met a documented acceptance criterion. The implant was released for transplantation based on the donor eligibility determination and a review of processing records.

Trace amounts of manufacturing residuals may remain after processing (acetic acid, acetone, hydrogen peroxide and sodium hydroxide).

STERILIZATION

Tutoplast is a tissue sterilization process that includes meticulous cleaning and gentle solvent dehydration of tissue. The process inactivates or removes potential pathogens, gently removes unwanted materials, such as cells, antigens and viruses and allows the implant to be stored at room temperature.

TUTOPLAST®



Low dose gamma irradiation is applied terminally to the product to achieve a sterility assurance level (SAL) of 10⁻⁶.

STORAGE AND SHIPPING



This symbol on the outer label indicates the storage temperature range for the product.



This symbol on the outer label indicates the expiration date of the product.

Store in a clean, dry environment at the temperature range specified on the product label. Keep away from sunlight.

SHIPPING CONDITIONS

Implant is shipped at ambient temperature via expedited shipping methods.

WARNINGS

The same potential medical/ surgical conditions or complications that apply to any surgical procedure may occur during or following implantation. The surgeon is responsible for informing the patient of the risks associated with their treatment and the possibility of complications or adverse reactions. As with any human tissue implant, it is not possible to guarantee freedom from transmission of infectious agents or other adverse reactions such as hypersensitivity, allergic or immune response.

PRECAUTIONS

Prior to use, the surgeon must become familiar with the implant and the surgical procedure. Poor general medical condition or any pathology that would limit the blood supply and compromise healing should be considered when selecting patients for procedures using this implant; as such conditions may compromise outcomes.

The implant should be used with caution in surgical sites where an active infection is present or in sites with poor perfusion. If the surgeon determines that the clinical circumstances require implantation in a site that is contaminated or infected, appropriate local and/or systemic anti-infective measures should be taken.

Appropriate placement and fixation of the implant are critical to success of the surgical procedure.

INSTRUCTIONS FOR USE

It is important to read and understand the following instructions prior to clinical use. Improper preparation technique may adversely affect the success of the surgical procedure.

GENERAL INSTRUCTIONS:

- Use on a single occasion for a single patient only. Once the package is opened, the implant must be used for the current procedure or discarded.
- The outermost packaging is non-sterile and is used to protect the implant during shipping and storage.
- Remove the double-barrier packaged product, the package insert, implant identification labels and Tissue Utilization Record from the box.
- Inspect the product, including all packaging and labeling materials carefully:
 - Do not use past expiration date specified on product label.
 - Do not use if the implant or packaging is damaged.
 - Do not use if there are discrepancies in label information.
- The implant's sterile barrier is comprised of two sealed pouches. To prevent contamination of the implant, use sterile technique for preparation and implantation.
- Additional product should be available in case of unexpected need during the procedure.
- Do not re-sterilize the implant.
- The implant and all packaging materials used by RTI Biologics, Inc. are latex-free.
- Use standard practices for handling and disposal of human tissue.
- Promptly report all product defects and patient adverse events to Coloplast (See Complaints and Returns section).

DIRECTIONS FOR IMPLANT PREPARATION:

1. Open outer pouch and pass inner pouch to sterile field.
2. Open sterile inner pouch and remove implant.
3. Re-hydrate the implant prior to use by soaking in sterile, room temperature saline solution for up to 30 minutes to improve suppleness and handling properties.
4. Pharmaceutical antibiotics or other antimicrobial agents prescribed by the surgeon as a precaution against incidental infection may be added to the soaking solution. The prescribing surgeon is responsible for selecting a suitable antibiotic or other antimicrobial agent at the appropriate concentration.
5. Size the implant according to the tissue defect. If necessary, multiple pieces of implant may be sutured together to cover a larger defect.
6. Place the implant securely to prevent displacement and to aid incorporation. Either absorbable or non-absorbable suture material may be used. Select the appropriate suture size for the surgical procedure. If absorbable sutures are used, it is recommended to select the longest lasting materials available. Place the stitches 2-3 mm from the edge of the implant. Use the implant where it is under minor to moderate tension.

TISSUE UTILIZATION RECORD (TUR CARD)

Complete and return the enclosed Tissue Utilization Record (TUR) to RTI Biologics, Inc. This information is kept confidential and used only for implant tracking. The TUR card should be submitted even if the implant was discarded. If information is not hand-written please verify labels or stamps are placed on both pages prior to mailing or faxing card. Refer to the enclosed TUR card for additional information.












CUSTOMER RETURNS AND COMPLAINTS

For further information or to report a complaint or adverse event, please contact:

Coloplast

1601 West River Rd. N
Minneapolis | MN 55411 | USA
Tel: 1.800.258.3476 | Fax: 1.866.216.4161
Website: www.us.coloplast.com

Axis is registered trademark of Coloplast.

Definition of Label Symbols		
		
Consult instructions for use	Expiration date	Storage temperature limits
		
Sterile by Gamma Irradiation	Single use only. Do not reuse	For prescription use only
		
Catalog number	Serial number (Tissue ID)	Lot number (Donor ID)
		
Keep Away from Sunlight	Manufacturer	

Manufactured By:

RTI Biologics, Inc.
11621 Research Circle | Alachua, FL 32615 | USA
Tel: 800.624.7238 or 386.418.8888
Fax: 386.462.5533
Email: cs@rtix.com

Tutoplast is a registered trademark of Tutogen Medical GmbH.