Digitex

Suture Delivery System

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CA Coloplast Canada Corporation
3300 Ridgeway Drive, Unit 12
Mississauga, Ontario L5L 5Z9
Canada
Tel: 1-905-820-7588
Fax: 1-905-820-8218

US Coloplast Corp.
1601 West River Road North
Minneapolis, MN 55411 USA
Toll Free: 800-328-3863
866-216-4161 fax
Direct: +1 612-337-7800
+1 612-337-7803 fax

70002133 Rev. C  **Coloplast** and **Digitex** are registered trademarks of Coloplast A/S. **Flexishaft** and **SutureSmart** are trademarks of Coloplast A/S. © 2012-03. All rights reserved Coloplast A/S, Humlebæk, Denmark.

Digitex® Suture Delivery System

SINGLE USE ONLY

Contents supplied sterile using Ethylene Oxide (EtO). Do not use if sterile barrier is damaged or appears damaged. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. If damaged, please contact your Coloplast representative.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Device Description

The **Digitex** Suture Delivery System (SDS) is composed of a delivery device and a suture cartridge and is designed for use by the physician to facilitate the placement of suture.

As demonstrated through data collected from a cadaveric study, the **Digitex** SDS can be used when direct visualization is not possible (i.e., placement is determined by palpation only), and/or when the anatomical location is difficult to reach. In addition, the device shaft incorporates a malleable portion, called **Flexishaft**, that has sufficient rigidity for functionality, but allows the device head to be repositioned 15 degrees in any direction. If desired, the user may adjust the device head position to more accurately target the desired anatomical landmark.

The **Digitex** Suture Cartridge is designed for use with the **Digitex** Delivery Device. **Digitex** Suture Cartridges are available in a variety of suture types and sizes. Detailed information regarding individual suture cartridge configurations can be found in the **Digitex** Suture Cartridge packaging. Each cartridge is preloaded with precut lengths of suture and are provided with either 0-0 or 2-0 suture. The suture has a specially designed cap at one end and is supplied with or without a crescent needle at the other end.

Indications

The **Digitex** Delivery Device is a sterile, disposable device intended to deliver a suture to the operative site. The device assists in suturing by passing a needle through the tissue to capture the suture, and suture ligation and knot-tying by holding the suture.

Contraindications

The use of The **Digitex** System is contraindicated in patients who are not candidates for surgical procedures and in applications requiring placement of suture into bone.

The **Digitex** Suture Cartridge is contraindicated in patients with known sensitivities or allergies to its components.

Warnings: Digitex Suture Delivery System

- The **Digitex** Suture Delivery System is intended for use only by physicians with adequate training and experience.
- The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures.

Warnings: Digitex Suture Cartridge

- Do not resterilize. Discard open, unused sutures.
- As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.
- Users should be familiar with surgical procedures and techniques involving suture before employing for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
- Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Precautions: Digitex Suture Delivery System

- As anatomy of individual patients may vary greatly, it is important that the intended locations for suture placement are planned for each procedure and each individual patient. Use of imaging methods may aid in confirming that non-targeted anatomical structures are not inadvertently injured.
- Use only the **Digitex** Suture Cartridge with the **Digitex** Delivery Device; the Delivery Device is designed to be used only with the **Digitex** Suture Cartridge which incorporates a specially designed suture cap.
- Do not activate the device without the cartridge and suture cap securely loaded. This may damage the operating mechanism for grasping the suture.
- After loading the cartridge onto the device, ensure that the cartridge is securely placed on its post and that the suture cap is contained within its receptacle.
- The **Digitex** SDS is not designed for placing suture into or through bone.

- The malleable shaft, **Flexishaft**, of the **Digitex** Delivery Device has been designed to function at various angles(+/- 15°), and the user should be aware of the resulting changes in needle delivery angle. Note: Subtle angular changes to the malleable **Flexishaft** are intended; changes greater than 15° may result in premature shaft fatigue.
- After confirming correct positioning, rock the Delivery Device head downward onto the tissue, applying light finger pressure to the base of device head. The device must be held firmly in position to drive the needle through the tissue and into the suture cap.
- The lever of the **Digitex** Delivery Device must be squeezed and completely actuated to properly seat the needle into the suture cap. Inadequate force at any of the fixation points could result in a failure to engage the needle into the cap.
- The **Digitex** SDS is not intended for use as a retractor, dissector, probe or purpose other than delivering a suture.
- Once the needle has been actuated, the **Digitex** Delivery Device should not be rotated, angled or otherwise moved to disrupt engagement of the needle and suture cap.
- If reposition of a suture is desired, clip off suture cap and needle, if applicable, and gently remove suture. Discard used needles appropriately in “Sharps” containers

Precautions: Digitex Suture Cartridge

In handling this or any other surgical suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Adequate knot security requires the accepted surgical technique of flat, squared ties as warranted by surgical circumstance and experience of the surgeon.

Adverse Events: Digitex Suture Delivery System

As with any instrument used for surgery, the potential complications associated with the use of the **Digitex** SDS to place sutures include:

- Bleeding
- Hematoma
- Infection
- Injury to vessels, nerves, internal organs and / or tissue including perforation and occlusion
- Inflammatory reaction to suture material or trauma of being sutured
- Inflammatory reaction to surgical procedure
- Dehiscence
- Erosion
- Extrusion
- Injury to the abdominal wall
- Peritonitis

Adverse Events: Digitex Suture Cartridge

Adverse effects associated with the use of the suture could include: wound dehiscence, calculous formation in the urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.

How Supplied: Digitex Delivery Device


The **Digitex** Delivery Device is packaged on a polypropylene support card in a double-barrier Tyvek/PET pouch. Four (4) packaged **Digitex** Delivery Devices are packaged in a shelf box. The devices will be shipped sterile (EtO sterilization).

How Supplied: Digitex Suture Cartridge

The **SutureSmart** cartridges are supplied sterile in precut lengths and prewound on a cartridge with and without an attached crescent needle. The suture is provided in either sizes 0-0 and 2-0 (USP). The suture cartridge is packaged individually in a single barrier Tyvek/PET pouch. Twelve pouches are packaged in a shelf box.

Instructions


1. Select desired suture size; cartridges are provided with either 0-0 or 2-0 suture.
2. Inspect **Digitex** Suture Delivery Device and **Digitex** Suture Cartridge to ensure components are not damaged. If any damage is noted, do not use the device and inform your Coloplast representative.
3. Delivery Device knob should be in “neutral position”; arrows should line up with wings in the horizontal position.
4. Suture Cartridge preparation: Grasp the suture cartridge. The clear side is the top (i.e., “up” side); grey is the bottom. Firmly seat cartridge on post. Unwind enough suture length to thread suture cap in device head. Insert suture cap into the device head, ensuring that the cap is fully seated. Wind any excess suture back onto cartridge by rotating the seated cartridge.

5. The Delivery Device shaft has a malleable portion that allows the device head to be repositioned 15° in any direction (approximately 5/8”). If repositioning of the head is desired to target the landmark, grasp the Delivery Device on each end of the malleable section and gently bend until the desired angle is achieved. The head of the device may then be positioned at the anatomical landmark to be sutured.


6. After confirming correct positioning, rock the head downward onto the tissue, applying light pressure to the base of the device head. Gently squeeze the lever, holding the Delivery Device firmly in position to drive the needle through the tissue and into the suture cap. Do not turn or twist the Delivery Device while activating. Care should be taken to apply little or no direct pressure on the distal curvature of the head during actuation.

7. A soft, tactile “click” will be felt when the needle is fully engaged with the suture cap. If a click is not felt, ensure that the lever has reached the full extent of its travel. Secure engagement may be possible even without the “click,” and this can be confirmed upon withdrawal of the Delivery Device.

8. Once engaged, ensure that the lever returns to its original, retracted position, and withdraw the Delivery Device from the field. The suture will unwind from cartridge and the suture cap will be retained within the Delivery Device upon withdrawal.

9. Once the Delivery Device head is well away from the tissue insertion point, rotate the wings on the release knob twice to disengage the suture cap from the Delivery Device needle.

10. Remove cartridge from Delivery Device and unwind suture so that cartridge may be set aside until it is time to tie the suture knots. Care should be taken when removing the full suture length from the cartridge as a crescent needle may be attached to the end of the suture. If desired, place additional sutures using the crescent

needle. If reposition of a suture is desired, clip off suture cap and crescent needle and gently remove suture. Do not pull suture cap back through tissue. Suture cap is not intended to be reused or implanted. Clip off suture cap and tie knots. Properly discard crescent needles in “Sharps” containers.

11. To place additional sutures, rotate the release knob twice to return the Delivery Device to actuation mode, rinse Delivery Device head in a sterile saline or antibiotic solution to remove particulates, and repeat steps 1-10, utilizing a new cartridge for each pass.

Product Evaluation

Coloplast requests physicians to notify the company of any complications which may develop with the use of this device, and requests return of used devices or components associated with the complication. For safe handling during shipment and upon receipt, Coloplast requests that devices be decontaminated prior to shipment. This is requested even though Coloplast will autoclave-sterilize any returned opened product. If necessary, Coloplast may analyze the device, and the patient and physician may be asked to allow Coloplast to perform tests that might alter the condition of the device.

Any complications from the use of this device should be brought to the immediate attention of your local sales representative's office. Contact information is located at the front of this booklet or by contacting: Quality Assurance, Product Evaluations Department, Coloplast Corp., 1601 West River Road North, Minneapolis, MN 55411. Toll-free telephone (800) 338-7908 in USA; or outside USA: (612) 337-7800

Return Goods Authorization

Authorization must be received from Coloplast prior to the return of merchandise. Product returned must have all manufacturer's seals intact and be returned within 30 days from date of invoice to be eligible for credit or replacement. Please contact the Coloplast Customer Service Department for details. To obtain a Return Authorization number contact: (800) 258-3476 or outside USA: (612) 337-7800; or fax (866) 216-4161 or outside USA: (612) 337-7803. Returned products may be subject to restocking charges.

Product Order Information

To order, please contact your local sales representative or Coloplast Customer Service Department at Coloplast, 1601 West River Road North, Minneapolis, MN 55411; Toll-free telephone: (800) 258-3476; or outside USA: (612) 337-7800; or fax (866) 216-4161 or outside USA: (612) 337-7803.

References

Literature references are available upon request from your local sales representative or: Coloplast Customer Service, 1601 West River Road North, Minneapolis, MN 55411 USA
Toll-free telephone: (800) 258-3476; or outside USA: (612) 337-7800

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