Restorelle® Y

Polypropylene mesh

Treillis en polypropylène

Restorelle® Y Polypropylene mesh

CAUTION: Federal (USA) law restricts this device to sale by, or on the order of a physician trained in performing synthetic mesh procedures for surgical repair of pelvic organ prolapse.

Restorelle® Y should only be used by surgeons familiar with the surgical procedures and techniques involving nonabsorbable meshes and who have adequate education and experience in the treatment of pelvic organ prolapse.

Your Coloplast sales representative can assist in identifying supplemental educational programs.

CAUTION: Read all information contained in this product label including Indications, Contraindications, Warnings, Precautions, Adverse Events, Instructions for Use, etc., prior to using this product.

Device Description:

Polypropylene mesh, using **Smartmesh®** and **Smartpores™** Technology.

Restorelle Y is a preshaped Y graft constructed of knitted non-absorbable monofilaments of polypropylene, a synthetic polymer.

Indications:

Restorelle Y is indicated for use as a bridging material for sacrocolposuspension/sacrocolpopexy (i.e., abdominal placement via laparotomy, laparoscopic, or robotic approach) where surgical treatment for vaginal vault prolapse is warranted.

Contraindications:

It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the contraindications associated with the use of this product. The product is contraindicated for use in the following conditions:

- Pregnancy or desire for future pregnancy
- Potential for further growth (e.g. adolescents)
- Pre-existing local or systemic infection. Treat the infection with the appropriate antiseptics and/or antibiotics to eliminate the infection before placing the **Restorelle Y** mesh.
- Any pathology, including known or suspected uterine pathology, which would compromise implant or implant placement.
- Taking anticoagulant therapy
- Sensitivity/allergy to polypropylene

Warnings and Precautions

It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the warnings and precautions associated with the use of this product and the associated surgical risks.

Warnings:

- **Restorelle Y** mesh should only be used by surgeons familiar with the surgical procedures and techniques involving non-absorbable meshes and who have adequate education and experience in the treatment of pelvic organ prolapse.
- Based on physician experience and education, a thorough assessment of each patient should be made to determine the suitability of a synthetic mesh procedure.
- The patient should be counseled that alternative non-mesh prolapse surgeries may be appropriate, and the reason for choosing a mesh procedure should be explained.
- Obtain patient consent prior to surgery and ensure that the patient has an understanding of the postoperative risks and potential complications of mesh surgery.
- Patient counseling should include a discussion that the mesh to be implanted is permanent, and that some complications associated with implanted mesh may require additional surgery; repeat surgery may not resolve these complications. Serious adverse tissue responses or infection may require removal of mesh.

Patient Warnings

Proper patient selection should be carefully considered in the following conditions as these factors may increase the likelihood of complications and adverse events:

- Age
- Smoking
- Connective tissue disorder
- Physical characteristics, (e.g. increased body mass index)
- Renal insufficiency or upper urinary tract obstruction
- Pelvic radiation therapy or chemotherapy in the last 12 months
- Uncontrolled diabetes, due to the potential risk of altered tissue healing

The risks and benefits of using **Restorelle** in patients with compromised immune systems and/or pathology that would limit blood supply and/or compromise healing should be carefully considered. The patient should be informed that any future pregnancy may negate the benefits of this surgical procedure. Patients should be instructed to immediately report any onset of bleeding, pain, abnormal vaginal discharge or signs of infection occurring at any time.

Product Warnings

- Do not use product that has damaged or opened packaging, as sterility may be compromised.
- Do not use product that has expired. Use of expired product can lead to infection or degradation and the product may not function as intended.
- **Restorelle Y** mesh is sold sterile for single use only, and should never be re-sterilized. In the event the product becomes contaminated prior to use, the device should be returned to Coloplast for replacement.
- Each device should be carefully examined prior to surgery and continuously monitored throughout the surgical procedure to ensure the structural integrity and sterility of the device has not been compromised in any way.
- A mesh which has been damaged or on which repairs have been attempted should not be implanted.
- If the mesh is removed from the patient, it should be handled and discarded as biologically hazardous.

Procedure Warnings

- The procedure to insert the **Restorelle Y** mesh requires good knowledge of pelvic anatomy.
- Good surgical practice must be followed during the **Restorelle Y** procedure.
- Prolapse may cause anatomical deformation.
- Avoid visceral perforation.
- Cystoscopy should be performed to confirm bladder and urethral integrity.
- A digital rectal exam should be performed to detect possible rectal perforation.
- Avoid placing excessive tension on the **Restorelle** mesh implant during placement and adjustment to maintain mesh integrity.
- Use of mesh may make future hysterectomy more difficult due to tissue in-growth and scarring.
- Additional consideration should be given to the use of **Restorelle Y** in any patient with a history of prior abdominal or pelvic surgeries. Consideration should also be given to the ability of the patient to tolerate the surgical procedure.
- There should be an appropriate margin of mesh extending beyond the fixation points.
- Inadequate suturing of the graft material to the pelvic tissue may lead to failure of the repair and recurrence of the prolapse.

Precautions:

Product

- Reuse of this single use product may create a potential risk to the patient. Reprocessing, cleaning, disinfection and/or resterilization may compromise a product characteristic which in turn creates an additional risk of physical harm to and/or infection of the patient.

Procedure

- Use caution to prevent intra-operative injury.
- Take care to avoid neurovascular injury. Observe patient for any signs of abnormal bleeding or clinical signs of nerve damage.
- Do not let **Restorelle Y** come into contact with sharp objects which could cause damage to the mesh (e.g., staples, clips, or clamps).

Postoperative Care:

Physician

- Standard post-operative protocols should be followed.
- Proper surgical practice should be followed for post-operative management of contaminated or infected wounds.
- General antibiotics should be administered at the physician's discretion.
- A catheter and/or vaginal pack should be used at the discretion of the physician.

Patient

- Physical strain, sexual intercourse and heavy lifting should be avoided for six weeks after surgery, but the patient may resume normal activities after two weeks or at the physicians discretion.
- Patients should be instructed to immediately report any onset of bleeding, pain, vaginal discharge or sign of infection occurring at any time.
- If infection occurs, partial or full mesh removal or revision may be necessary, per surgeon discretion.

CA Coloplast Canada Corporation
3300 Ridgeway Drive, Unit 12
Mississauga, Ontario L5L 5Z9
Canada
Tel: 1-905-820-7588
Fax: 1-905-820-8218

US Coloplast Corp.
1601 West River Road North
Minneapolis, MN 55411 USA
Toll Free: 800-328-3863
866-216-4161 fax
Direct: +1 612-337-7800
+1 612-337-7803 fax

Coloplast

Manufacturer / Fabricant :

Coloplast A/S
Holtedam 1
3050 Humlebaek
Denmark

Adverse Events:

Potential adverse events are those associated with surgery using implantable synthetic mesh materials. As with all foreign bodies, **Restorelle** mesh is likely to exacerbate any existing infection. Local irritation at the wound site and/or a foreign body response may occur. There is also the risk of complete failure of the procedure resulting in recurrent prolapse. The following complications are known to occur with synthetic mesh implantation:

- Mesh erosion (e.g., vaginal, urethral, bladder)
- Mesh extrusion
- Mesh exposure
- Infection
- Pain (acute or chronic)
- Bladder, bowel, urethra, vagina, vessel, and/or nerve perforation/injury

Additional known risks in prolapse surgical procedures with synthetic mesh may include but are not limited to:

- Allergic and/or foreign body reaction (e.g., granuloma formation)
- Adhesion and/or abscess formation
- Bleeding
- Chronic inflammation and/or discharge
- Defecatory dysfunction
- Delayed wound healing
- De novo and/or worsening dyspareunia
- De novo, continued, and/or worsening incontinence
- De novo or worsening prolapse in untreated compartment
- De novo overactive bladder
- Dyspareunia/partner discomfort
- Fecal incontinence and/or anal sphincter incompetence
- Fistula formation
- Hematoma, seroma, and/or hemorrhage
- Ileus
- Need for blood transfusion
- Neurologic and/or neuromuscular symptoms
- Palpable mesh (patient and/or partner)
- Sexual dysfunction
- Small bowel obstruction
- Suture erosion
- Transient or permanent urinary retention/obstruction
- Urinary tract infection
- Vaginal bulge
- Vaginal rigidity
- Vaginal scarring and/or tissue contracture
- Voiding dysfunction
- Wound dehiscence and/or necrosis
- Wound /surgical site infection

The occurrence of these events may require partial or complete removal of the mesh. Patients should be monitored regularly after the device has been implanted for immediate treatment of any adverse event.

Instructions for Use:

The mesh can be used as an inlay to support fascial repair. Sutures or other means may be used to fix the mesh in place. Specific procedures may vary according to judgments of skilled medical practitioners and variations in the pelvic characteristics of individual patients. Carefully examine the device prior to surgery and continuously monitor throughout the surgical procedure to ensure the structural integrity of the device is not compromised in any way. Good surgical practice must be followed during the **Restorelle** procedure.

1. Presacral Space Dissection. Expose anterior longitudinal ligament for later graft fixation
2. Dissect the Vesicovaginal and Rectovaginal Spaces. Expose vaginal apex for graft suture fixation
3. Vaginal Apex Graft Attachment. Size the mesh to the individual patient anatomy
4. Anterior Longitudinal Ligament Graft Attachment. Size the mesh to the individual patient anatomy
5. Peritoneal Closure. Retroperitonealize the graft
 - Inadequate suturing of the graft material to the pelvic tissue may lead to failure of the repair and recurrence of the prolapse.
 - There should be an appropriate margin of mesh extending beyond a suture line. Inadequate fixation of the graft material to the pelvic tissue may lead to failure of the repair and recurrence of the prolapse.
 - Avoid excessive tension on the mesh implant during placement.
 - It is recommended to confirm bladder and rectal integrity with cystoscopy and digital rectal exam.

Product Complaints:

Coloplast requests physicians to notify the company of any complications which may develop with the use of this device, and requests return of any explanted devices or components. For safe handling during shipment and upon receipt, Coloplast requests that devices be decontaminated prior to shipment. This is requested even though Coloplast will autoclave-sterilize any opened product returned. Alteration for the purposes of venting to prevent additional damage will be performed as required. If explant is necessary, Coloplast may analyze the explanted device, and the patient and physician may be asked to allow Coloplast to perform tests that might alter the condition of the device.

Any complications from the use of this device should be brought to the immediate attention of your local sales representative's office. Contact information is located in this booklet or by contacting: Quality Assurance, Product Evaluations Department, Coloplast, 1601 West River Road North, Minneapolis, MN 55411. Toll-free telephone: (800) 338-7908 in USA; or outside USA: (612) 337-7800.

Returned Goods Authorization:

Authorization must be received from Coloplast prior to the return of merchandise. Merchandise returned must have all manufacturer's seals intact and be returned within 30 days from date of invoice to be eligible for credit or replacement. Please contact the Coloplast Customer Service Department for details. To obtain a Return Authorization number, see contact information located in this booklet or by contacting: (800) 258-3476 or outside USA: (612) 337-7800; or fax (866) 216-4161 or outside USA: (612) 337-7803. Returned products may be subject to restocking charges.

Product Order Information:

To order, please contact your local sales representative or Coloplast Customer Service Department at Coloplast, 1601 West River Road North, Minneapolis, MN 55411. Toll-free telephone: (800) 258-3476; or outside USA: (612) 337-7800; or fax (866) 216-4161 or outside USA: (612) 337-7803.

References:

Literature references are available upon request from your local sales representative or: Coloplast Customer Service, 1601 West River Road North, Minneapolis, MN 55411 USA. Toll-free telephone: (800) 258-3476; or outside USA: (612) 337-7800.

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