Supris®

Retropubic Sling System

CA Coloplast Canada Corp.
3300 Ridgeway Drive, Unit 12
Mississauga, Ontario L5L 5Z9
Canada
Tel: 1-905-820-7588
Fax: 1-905-820-8218

US Coloplast Corp.
1601 West River Road North
Minneapolis, MN 55411 USA
Toll Free: 800-328-3863
866-216-4161 fax
Direct: +1 612-337-7800
+1 612-337-7803 fax

Supris® Retropubic Sling System

CAUTION: Federal (U.S.A.) law restricts this device to sale by, or on the order of a physician trained in performing synthetic sling procedures.

Warning: Cystoscopy should be performed to confirm bladder integrity or to recognize bladder perforation.

Contact your local sales representative to arrange for educational opportunities with the Supris Retropubic Sling system. **Device Description:** The Coloplast **Supris** sling is a permanent, synthetic sub-urethral sling that is provided with disposable introducer needles in the **Supris** Retropubic Sling System. It is indicated for the surgical treatment of all types of stress urinary incontinence (SUI), for female stress urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency (ISD). **Supris** sling is made from knitted, monofilament polypropylene and has low elasticity. This structure gives the **Supris** sling resistance to traction, allows tissue colonization and facilitates positioning during surgery.

The **Supris** sling is provided sterile and for single-use only. The **Supris** sling is sterilized by ethylene oxide. Sterile, single-use only, introducers are included in this system. The introducers are manufactured from polypropylene (handles) and medical grade stainless steel (Introducer Needles) and are sterilized by ethylene oxide. These introducers can only be used with Coloplast **Supris** sling.

Indications: The **Supris** Retropubic Sling System is an implantable, sub-urethral, support sling indicated for the surgical treatment of all types of stress urinary incontinence (SUI), for female urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency. The sling is placed retropubically using two disposable introducers using either a “top down” or “bottom up” surgical approach.

Contraindications: It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the contraindications associated with the use of this product. The product is contraindicated for use in the following women:

- Pregnancy or desire for future pregnancy
- Potential for further growth (e.g. adolescents)
- Known active urinary tract infection and/or infection in operative field
- Sensitivity/allergy to polypropylene
- Taking anti-coagulant therapy

Warnings and Precautions

It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the possible warnings associated with the use of this product and the associated surgical risks.

Warnings:

- Do not use product that has damaged or opened packaging, as sterility may be compromised.
- The **Supris** Retropubic Sling System is sold sterile for single-use only, and should never be resterilized. In the event the product becomes contaminated prior to use, the device should be returned to Coloplast for replacement.
- Each device should be carefully examined prior to surgery and continuously monitored throughout the surgical procedure to ensure the structural integrity and sterility of the device has not been compromised in any way. A device which has been damaged or on which repairs have been attempted should not be used or implanted.
- The risks and benefits of using **Supris** in patients with compromised immune systems or any other conditions that affect healing should be carefully considered.
- The **Supris** Retropubic Sling System should only be used by surgeons who are qualified to perform this type of surgery and who are familiar with the use of non-absorbable mesh and the specific insertion technique for **Supris**.
- The procedure to insert the device requires a good knowledge of local anatomy and the correct use of the Introducer Needles in order to avoid damage to adjacent anatomical structures.
- An introducer that has been damaged should not be repaired or utilized.
- Avoid perforation of the viscus. Cystoscopy should be performed after placement to confirm bladder and urethral integrity.
- Avoid placing excessive tension on the **Supris** sling during placement and adjustment to maintain sling integrity and to avoid compression of the urethra when tensioning.
- Avoid placing **Supris** in patients with an abnormal urethra (e.g. fistula, diverticulum).
- The patient should be informed that any future vaginal delivery will negate the benefits of this surgical procedure and the patient is likely to become incontinent again.
- Consideration should be given to patients with uncontrolled diabetes, for the potential risk of altered tissue healing.
- If the product is removed, it should be handled and discarded as biologically hazardous.

Precautions:

- Reuse of this single-use product may create a potential risk to the user. Reprocessing, washing, disinfection and/or sterilization may compromise a product characteristic, which in turn creates an additional risk of physical harm to or infection of the user.
- Good surgical practice must be followed during the **Supris** Retropubic Sling System procedure.
- Special care should be taken in cases of bladder prolapse because of anatomical deformation. If the patient requires a cystocele repair, it should be completed prior to using the **Supris** sling.
- Patients who smoke should be carefully considered.
- Patients with connective tissue disorder should be carefully considered.
- Physicians should consider physical characteristics (e.g. body mass index) in determining who is a candidate for the **Supris** procedure.
- Use caution to prevent intraoperative injury.
- Use with caution in patients with renal insufficiency or upper urinary tract obstruction.
- Use with caution in any patient who has received pelvic radiation therapy or who has received chemotherapy in the last 12 months.
- Take care to avoid nonvascular perforation. Observe patient for any signs of retropubic/periurethral bleeding or clinical signs of nerve damage inter/postoperatively.
- Do not let the **Supris** sling come into contact with sharp objects that could cause damage to the device (i.e., staples, clips, or clamping on the sling).
- Patients should immediately report any onset of bleeding, pain, vaginal discharge or sign of infection that occurs at any time.

Adverse Events: Potential adverse events are those associated with surgery using implantable synthetic mesh materials. As with all foreign bodies, the **Supris** sling is likely to exacerbate any existing infection.

Transitory local irritation at the wound site and/or a foreign body response may occur. There is also the risk of complete failure of the procedure resulting in continued incontinence due to incomplete support or

overactive bladder. Known risks of incontinence procedures with synthetic slings may include, but are not limited to:

- Allergic and/or foreign body reaction (e.g. granuloma)
- Chronic inflammation and/or discharge
- De novo and/or worsening dyspareunia
- De novo overactive bladder
- De novo urge incontinence
- Dyspareunia/partner discomfort
- Erosion (e.g. vaginal, urethral)
- Fistula formation
- Hematoma and/or hemorrhage
- Impact on sexual function
- Infection
- Mixed urinary incontinence
- Need for a blood transfusion
- Nerve injury
- Neurologic symptoms
- Pain (acute or chronic)
- Sling migration
- Transient or permanent urinary retention/obstruction
- Urethral obstruction
- Urinary tract infection
- Vaginal contracture
- Vaginal exposure and/or extrusion
- Vaginal scarring
- Vascular injury
- Viscus perforation/injury
- Voiding dysfunction
- Worsening incontinence
- Wound dehiscence and/or necrosis
- Wound/surgical site infection

The occurrence of these events may require partial or complete removal of the sling. Patients should be monitored regularly after the device has been implanted for immediate treatment of any adverse event.

Instructions for Use:

Supris Retropubic (Top Down) Surgical Procedure
Surgery Preparation

1. Conduct proper anatomical evaluation of the patient.
2. Apply appropriate local, regional or general anaesthesia to the patient.
3. Insert Foley catheter in urethra and empty the bladder.
4. Place patient in lithotomy position.
5. Use suitable vaginal retraction.

Surgical Steps

1. Make a full thickness vertical incision on the anterior vaginal wall at the junction of the middle and lower third of the urethra.
2. Perform blunt paraurethral dissection. Allow enough room for your finger to meet the tip of the needle and manually guide the introducer through the vaginal incision.
3. Make two abdominal skin incisions approximately 1 cm in length and 2 cm from the midline or one finger breadth above and lateral to the pubic symphysis. Avoid incisions that may go too far laterally.
4. Insert the needle into the abdominal incisions. Carefully press the needle down until it passes through the anterior rectus fascia and passes along the back of the symphysis pubis. Rotate the needle vertically to the posterior side of the pubic bone.
5. Move the needle along the posterior side of the pubic bone toward the vaginal incision. Use your finger to guide the needle through the endopelvic fascia through the vaginal incision.
6. Repeat step 5 on the opposite side with the second introducer needle.
7. Cystoscopy should be performed after the passing of the needles to confirm bladder integrity or to recognize bladder perforation.
8. Attach **Supris** sling to the introducers by threading the sling through the eye of the introducer. Pull the sling through the eyelet about 3-4 centimeters. Reverse the introducer through the incision tunnel to place the sling.
9. Vaginal retraction should be removed prior to tensioning. To prevent over-tensioning, there should be visible space between the urethra and sling. Make sure that the sling is lying flat under the urethra.
10. Cut off excess sling at the abdominal incisions. Push the skin down without pulling on the sling. Cut the end of the **Supris** sling below the level of the skin. The sling should lie within the subcutaneous tissue.
11. Carefully close the vaginal incision. Close the abdominal incisions.

Supris Retropubic (Bottom Up) Surgical Procedure
Surgery Preparation

1. Conduct proper anatomical evaluation of the patient.
2. Apply appropriate local, regional or general anaesthesia to the patient.
3. Insert Foley catheter in urethra and empty the bladder.
4. Place patient in lithotomy position.
5. Use suitable vaginal retraction.

Surgical Steps

1. Make a full thickness vertical incision on the anterior vaginal wall at the junction of the middle and lower third of the urethra.
2. Perform blunt paraurethral dissection. A paraurethral tunnel is developed up to the underside of the endopelvic fascia on either side of the urethra by scissors and/or blunt dissection.
3. Make two abdominal skin incisions approximately 1 cm in length and 2 cm from the midline or one finger breadth above and lateral to the pubic symphysis. Avoid incisions that may go too laterally.
4. Attach **Supris** sling to the introducers by threading the sling through the eye of the introducer. Pull the sling through the eyelet about 3-4 centimeters.
5. Insert the introducer into the vaginal incision. The introducer is passed in the paraurethral tunnel, staying close to the inferior pubic ramus, penetrating the endopelvic fascia, passing through the retropubic space and then the rectus complex to exit the skin just above the pubic tubercle 2 cm from the midline. Care should be taken to avoid bladder perforation.

- Cystoscopy should be performed after the passing of the introducer to confirm bladder integrity or to recognize bladder perforation.
- Reverse the introducer through the incision tunnel to place the sling.
- Repeat passage on the opposite side with the second introducer.
- Cystoscopy should again be performed after the passing of the second introducer to confirm bladder integrity or to recognize bladder perforation.
- Reverse the introducer through the incision tunnel to place the sling.
- Vaginal retraction should be removed prior to tensioning. To prevent over-tensioning, there should be visible space between the urethra and sling. Make sure that the sling is lying flat under the urethra.
- Cut off excess sling at the abdominal incisions. Push the skin down without pulling on the sling. Cut the end of the **Supris** sling below the level of the skin. The sling should lie within the subcutaneous tissue.
- Carefully close the vaginal incision. Close the abdominal incisions.

Postoperative Care:

- Standard postoperative protocols should be followed, including standard practice for postoperative management of contaminated or infected wounds.
- General antibiotics should be administered at the physician's discretion.
- A catheter and/or vaginal pack should be used at the discretion of the surgeon.
- Physical strain, sexual intercourse and heavy lifting should be avoided for six weeks after surgery, but the patient can resume other normal activities after two weeks at the surgeon's discretion.
- Patients should immediately report any onset of bleeding, pain, vaginal discharge or sign of infection that occurs at any time.
- If infection occurs, partial or full sling removal or revision may be necessary, per physician discretion.

Product Complaints: Coloplast requests physicians to notify the company of any complications which may develop with the use of this device, and requests return of any explanted devices or components. For safe handling during shipment and upon receipt, Coloplast requests that devices be decontaminated prior to shipment. This is requested even though Coloplast will autoclave-sterilize any opened product returned. Alteration for the purposes of venting to prevent additional damage will be performed as required. If explantation is necessary, Coloplast may analyze the explanted device, and the patient and physician may be asked to allow Coloplast to perform tests that might alter the condition of the device.

Any complications from the use of this device should be brought to the immediate attention of your local sales representative's office. Contact information is located in this booklet or by contacting: Quality Assurance, Product Evaluations Department, Coloplast,1601 West River Road North, Minneapolis, MN 55411. Toll-free telephone: (800) 338-7908 in USA; or outside USA: (612) 337-7800.

Returned Goods Authorization: Authorization must be received from Coloplast prior to the return of merchandise. Merchandise returned must have all manufacturer's seals intact and be returned within 30 days from date of invoice to be eligible for credit or replacement. Please contact the Coloplast Customer Service Department for details. To obtain a Return Authorization number, see contact information located in this booklet or by contacting: (800) 258-3476 or outside USA: (612) 337-7800; or fax (866) 216-4161 or outside USA: (612) 337-7803. Returned products may be subject to restocking charges.

Product Order Information: To order, please contact your local sales representative or Coloplast Customer Service Department at Coloplast, 1601 West River Road North, Minneapolis, MN 55411. Toll-free telephone: (800) 258-3476; or outside USA: (612) 337-7800; or fax (866) 216-4161 or outside USA: (612) 337-7803.

References: Literature references are available upon request from your local sales representative or: Coloplast Customer Service, 1601 West River Road North, Minneapolis, MN 55411 USA. Toll-free telephone: (800) 258-3476; or outside USA: (612) 337-7800.

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