

CLEANING BRUSHES FOR SIALENDOSCOPES



Don't Trust Your Sialendoscopes with Anything Less

Designed for effective and efficient cleaning, our single-use brushes ensure compliance with product IFUs and help maintain the integrity and precision of our endoscopes.

- The only brushes designated for use with KARL STORZ sialendoscopes
- Gentle, thorough, confident cleaning
- Available in various shapes and sizes
- 10% discount on convenient bundled packages

STORZ
KARL STORZ – ENDOSKOPE

CLEANING BRUSHES FOR SIALENDOSCOPES






Don't Trust Your Sialendoscopes with Anything Less

- KARL STORZ single-use cleaning brushes are specifically designed for use with KARL STORZ endoscopes—making it easy to comply with our product IFUs
- The nonabrasive bristles provide a soft yet effective scrubbing action while protecting delicate endoscopes from damage
- The comprehensive range of styles and sizes covers our complete line of sialendoscopes

ORDERING INFORMATION FOR BRUSH BUNDLE PACKAGES

SIALBRUSH-KT1	SIALBRUSH-KT2	SIALBRUSH-KT3	SIALBRUSH-KT4
<ul style="list-style-type: none"> • 6 PACKS: 27652-05 • 3 PACKS: 27651K4 	<ul style="list-style-type: none"> • 6 PACKS: 27652-05 • 3 PACKS: 27651K1 • 3 PACKS: 27651K4 	<ul style="list-style-type: none"> • 6 PACKS: 27652-05 • 3 PACKS: 27651K2 • 3 PACKS: 27651K4 	<ul style="list-style-type: none"> • 6 PACKS: 27652-05 • 3 PACKS: 27651K3 • 3 PACKS: 27651K4
USE WITH: 0.89 (11581A) 0.8 (11572A)	USE WITH: 1.1 (11582A) 1.1 (11573A)	USE WITH: 1.3 (11575A) 1.6 (11574A)	USE WITH: 1.6 (11583A)



Featured Brushes			
	27652-05	Single use flat cleaning brush, for the exterior cleaning. 5 per package. 15 mm width x 15.1 cm length x 7.5 mm height.	All sialendoscopes
	27651K1	Single use cleaning brush, for working channels with a diameter of 0.4 – 0.6 mm. 10 per package. 40 cm in length.	1.1 (11582A) & 1.1 (11573A)
	27651K2	Single use cleaning brush, for working channels with a diameter of 0.6 – 0.8 mm. 10 per package. 40 cm in length.	1.3 (11575A) & 1.6 (11574A)
	27651K3	Single use cleaning brush, for working channels with a diameter of 0.8 – 1.4 mm. 10 per package. 40 cm in length.	1.6 (11583A)
	27651K4	Single use cleaning brush, for working channels with a diameter of 0.2 – 0.4 mm. 10 per package. 40 cm in length.	All sialendoscopes

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com