



Introducing BD Spirit™ Leg Bags & Drainage Bags

Customizable tubing length

Leg Bags include a freestanding connector, allowing customization of the tubing length for optimal placement.

BD® Spirit™ Leg Bag

Soft non-woven backing

For all-day wear.

Ergonomically designed

Comprised of three distinct chambers compared to traditional single-chamber designs.

Ridged connector

Compatible with BD® Male External Catheters and Foley catheters.

Leg straps

One pair of fabric leg straps featuring silicone wave lining included with every leg bag.

Tap-Valve

Preconnected to silicone tubing for connection to an overnight system.



Introducing BD Spirit™ Collection Systems

Collection devices:

Previous Bag Order Codes	Previous Bag Description	Current Bag Order Codes	Current Bag Description
150106	9 oz. Leg bag with rubber cap valve	USLB350S	350mL Spirit™ Leg Bag Tap-Valve
150719	19 oz. Leg bag with Flip-Flo™ valve	USLB500S	500mL Spirit™ Leg Bag Tap-Valve
150102	19 oz. Leg bag rubber cap valve		
150819	19 oz. Leg bag with Flip-Flo™ valve, 18" PVC extension tubing		
150732	32 oz. Leg bag with Flip-Flo™ valve	USLB750S	750mL Spirit™ Leg Bag Tap-Valve
150103	32 oz. Leg bag rubber cap valve		
150832	32 oz. Leg bag with Flip-Flo™ valve, 18" PVC extension tubing		
154002	Bard® Urinary Drainage Bag	USBB2000S	2000mL Spirit™ Bed Bag Tap-Valve
802001	Bardia™ Drainage Bag		

For more information contact your BD sales representative.



Indications for use:

Spirit™ Leg Bag: Leg bags are indicated for urine collection with urinary catheters as determined by a physician.

Spirit™ Night/Bed Drainage Bag: Spirit™ Night/Bed Drainage Bags are indicated for urine collection for patients requiring urinary catheterization.

Warnings:

Spirit™ Leg Bag: Remove leg straps for conditions including, but not limited to: discoloration, swelling, numbness, rash, discomfort or if the leg strap is soiled. For pre-existing conditions, and/or if experiencing problems with use of the device, please consult your healthcare professional for assistance.

Spirit™ Night/Bed Drainage Bag: This is a single use device. Do not resterilize any portion of this device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and with applicable laws and regulations.

Please consult package insert for detailed safety information and complete instructions for use

BD, Covington, GA, 30014, U.S.

bd.com

BD, the BD Logo, Bard, Bardia and Spirit are trademarks of Becton, Dickinson and Company or its affiliates.
© 2025 BD. All rights reserved. 2510-41 BD-164964

