

ENDOMAT® SELECT/EQUIMAT® FLUID-MANAGEMENT SYSTEM



The Choice is Clear

The ENDOMAT® SELECT pump, together with the EQUIMAT® controller, offers precise fluid management through a simple user interface. This system includes a pressure controlled infusion pump, a fluid collection system, and a fluid deficit monitoring system. The combination is ideal for users seeking:

- Maximum patient safety
- Adjustable flow and pressure rates for optimal uterine distention
- Simplified setup and operation
- Streamlined workflow for the OR

ENDOMAT® SELECT/EQUIMAT® FLUID-MANAGEMENT SYSTEM

The Choice is Clear

- ① The **EQUIMAT®** fluid-monitoring system precisely measures and monitors fluid deficit in real time for patient safety. The device features visual and audible status indicators and warns the user when selected fluid deficit limit values are exceeded. Using the “Reset” button, the irrigation bags can be easily exchanged, even during the surgery.
- ② The **ENDOMAT® SELECT** is a single-roller pump for the irrigation of fluids. The unit combines ease of use and a clearly arranged user interface for greater user comfort, maximum patient safety, and optimized workflow.

Ordering Information

ENDOMATSELECTFM-KT Kit Includes:

UP210-GYN	ENDOMAT® SELECT Pump
20302020-1	EQUIMAT® (SCB)
38332130	Scale Measuring Element
26600EQC	Calibration Test Weight (500 g)
26602EQ	Fluid-Management Stand Set (includes: 26602EQB, 26602EQP, 20302037)

Software

UP602	Hysteroscopy Software , license (allows selection of “HYS” procedure)
UP610	Advanced Paket-Software , license (extends functions of installed software packages)

Accessories

031523-10-GYN	Tubing Set , irrigation, pressure-controlled, sterile, for single use, package of 10
030247-10	Single-use Y-Tubing Set , for suction, sterile, for single use, package of 10
26163FD	Kimberly Clark Collection Drape , sterile, for single use, 1 pk



Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com