

# KARL STORZ C-MAC® Blade Container System



## Protect Your C-MAC® Blade Investment

Our KSZ-39301CMAC Container System is specially designed to protect two C-MAC® blades during transport, reprocessing and storage.

- Houses two C-MAC® video laryngoscopes
- Save time and money over traditional blue-wrap trays
- Anodized to prevent pitting and corrosion


**STORZ**  
KARL STORZ—ENDOSKOPE




# KARL STORZ C-MAC® BLADE CONTAINER SYSTEM

## Protect Your C-MAC® Blade Investment

The KSZ-39301CMAC blade sterilization container system is designed to house your KARL STORZ C-MAC® video laryngoscopes through transport, sterilization and storage. It offers hassle-free, cost-effective sterilization for your valuable C-MAC® blade investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They help eliminate the need for blue sterilization wrap, while saving time and money.

Validated and approved for STERRAD® and V-PRO®\* sterilization systems\*

Ordering Information		
	KSZ-39301CMAC	<b>Container System</b> for two reusable C-MAC® Video Laryngoscope Blades. Includes container base, basket, and lid. 12.3" x 11" x 4"

Accessories for Vaporized H <sub>2</sub> O <sub>2</sub> /STERRAD®/V-PRO® Sterilization		
	KSZ-LTRNDFLTR	<b>Disposable Polypropylene Filters</b> , round, SteriTite® polypro filter for vaporized H <sub>2</sub> O <sub>2</sub> /STERRAD®/V-PRO® sterilization. 500 per package. 7.5" diameter.
	KSZ-LTCRDLG	<b>Large Load Cards</b> , process indicator cards for containers, low temperature, SteriTite® load cards for vaporized H <sub>2</sub> O <sub>2</sub> /STERRAD®/V-PRO® sterilization. 500 per package
	KSZ-LTSEAL	<b>Tamper-Evident Seals</b> , low temperature, for vaporized H <sub>2</sub> O <sub>2</sub> /STERRAD®/V-PRO® sterilization. 500 per package.

\*STERRAD® is a registered trademark of Advanced Sterilization Products. V-PRO® is a registered trademark of STERIS Corp.

\*\* Please refer KSZ-IFU\_SteriTite IFU KARL STORZ for more details on sterilization modalities. Refer to the product IFU to determine correct sterilization modality prior to reprocessing.

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue | El Segundo, CA 90245-5017, USA  
Phone: 424-218-8100 Toll-free: 800-421-0837  
Fax: 424-218-8525 Toll-free: 800-321-1304  
E-Mail: [communications@karlstorz.com](mailto:communications@karlstorz.com)

[www.karlstorz.com](http://www.karlstorz.com)

FOLLOW KARL STORZ



Caution: Federal (USA) law restricts this device to sale by or on the order of physician.  
All trademarks and trade names are the property of their respective owners.