


KARL STORZ LOTTA® AND Little LOTTA® Telescope Container System



Protect Your LOTTA® and Little LOTTA® Telescope Investment

Our KSZ-39301LTA Sterilization Container System is specially designed to protect your LOTTA® telescopes and accessories during transport, reprocessing and storage.

- Validated and approved for steam sterilization
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion





STORZ
KARL STORZ — ENDOSKOPE

KARL STORZ LOTTA® AND Little LOTTA® Telescope Container System

Protect Your LOTTA® Telescope Investment

The KSZ-39301LTA Sterilization Container System is designed to house your KARL STORZ LOTTA® and Little LOTTA® telescope and accessories through transport, sterilization and storage. It offers cost-effective, hassle-free sterilization for your valuable LOTTA® investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They also help eliminate the need for blue sterilization wrap, while saving time and money.

Ordering Information			
		KSZ-39301LTA	Container System, for LOTTA® or Little LOTTA® scopes and accessories. Includes basket, container base, and lid. 17" x 11" x 4.5"
		KSZ-39316F	Container System, for LOTTA® instruments. Includes two instrument racks, base basket, container base, and lid. 23.5" x 11" x 8"

Accessories For Steam Sterilization		
	KSZ-STRNDFLTR	Disposable Paper Filters , round, for steam sterilization. 100% cellulose, water repellent. 500 per package. 7.5" diameter.
	KSZ-STCRDLG	Large Load Cards , dual process indicator cards for steam sterilization. ID card with a dual-indicator dot. 500 per package.
	KSZ-STSTRIP	Dual Chemical Indicator Strips , with room for processing identification, 8" long with perforation to reconfigure up to two 4" long indicator strips. 500 per package.
	KSZ-STSEAL	Tamper-Evident Seals , for steam sterilization. 500 per package.

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue | El Segundo, CA 90245-5017, USA
Phone: 424-218-8100 Toll-free: 800-421-0837
Fax: 424-218-8525 Toll-free: 800-321-1304
E-Mail: communications@karlstorz.com

www.karlstorz.com

FOR U.S. USE ONLY

FOLLOW KARL STORZ



Caution: Federal (USA) law restricts this device to sale by or on the order of physician.
All trademarks and trade names are the property of their respective owners.