

Cleaning Brushes For Sialendoscopes



Don't Trust Your Sialendoscopes with Anything Less

Designed for effective and efficient cleaning, our single-use brushes ensure compliance with product IFUs and help maintain the integrity and precision of our endoscopes.

- The only brushes designated for use with KARL STORZ sialendoscopes
- Gentle, thorough, confident cleaning
- Available in various shapes and sizes
- 10% discount on convenient bundled packages

STORZ
KARL STORZ — ENDOSKOPE







Cleaning Brushes For Sialendoscopes

Don't Trust Your Sialendoscopes with Anything Less

- KARL STORZ single-use cleaning brushes are specifically designed for use with KARL STORZ endoscopes—making it easy to comply with our product IFUs
- The nonabrasive bristles provide a soft yet effective scrubbing action while protecting delicate endoscopes from damage
- The comprehensive range of styles and sizes covers our complete line of sialendoscopes

All kits come with 30 of each brush.

SIALBRUSH-KT1	SIALBRUSH-KT2	SIALBRUSH-KT3	SIALBRUSH-KT4	SIALBRUSH-KT5
27652-05 27651K5	27652-05 27651K1 27651K4	27652-05 27651K2 27651K5	27652-05 27651K3 27651K5	27652-05 27651K1 27651K5
USE WITH: 0.89 (11581A) 0.8 (11572A)	USE WITH: 1.1 (11582A)	USE WITH: 1.3 (11575A) 1.6 (11574A)	USE WITH: 1.6 (11583A)	USE WITH: 1.1 (11573A)

Featured Brushes			
	27652-05	Single-Use Flat Cleaning Brush , for the exterior cleaning. 5 per package. 15 mm width x 15.1 cm length x 7.5 mm height	All sialendoscopes
	27651K1	Cleaning Brush , round, flexible, for working channel diameter 0.4–0.6 mm, length 40 cm, for single use, package of 10.	1.1 (11582A), 1.1 (11573A)
	27651K2	Cleaning Brush , round, flexible, for working channel diameter 0.6–0.8 mm, length 40 cm, for single use, package of 10.	1.3 (11575A), 1.6 (11574A)
	27651K3	Cleaning Brush , round, flexible, for working channel diameter 0.8–1.4 mm, length 40 cm, for single use, package of 10.	1.6 (11583A)
	27651K4	Cleaning Brush , round, flexible, for irrigation channel diameter 0.2–0.25 mm, length 30 cm, for single use, package of 10.	1.1 (11582A)
	27651K5	Cleaning Brush , round, flexible, for irrigation channel diameter 0.25–0.4 mm, length 40 cm, for single use, package of 10.	0.89 (11581A), 0.8 (11572A), 1.1 (11573A), 1.3 (11575A), 1.6 (11574A), 1.6 (11583A)

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue | El Segundo, CA 90245-5017, USA
Phone: 424-218-8100 Toll-free: 800-421-0837
Fax: 424-218-8525 Toll-free: 800-321-1304
E-Mail: communications@karlstorz.com

www.karlstorz.com

FOR U.S. USE ONLY

FOLLOW KARL STORZ



Caution: Federal (USA) law restricts this device to sale by or on the order of physician.
All trademarks and trade names are the property of their respective owners.