



Expanding the Boundaries of **ENDOLUMINAL THERAPIES**



DiLumen[®] EZ¹ **Endoluminal Interventional Platform**

Building on the established innovations of the DiLumen Platform, the DiLumen EZ¹ provides a reliable solution to support advanced endoluminal procedures in the lower gastrointestinal tract.

The single-balloon configuration offers unmatched stability, procedural precision, visualization, and safety — even under anatomically challenging conditions.

The overtube, featuring a hydrophilic coating, provides flexible guidance with superior gliding characteristics.

Providing essential benefits for endoluminal procedures



Confident navigation

- Unique one-balloon design for improved stability
- Atraumatic conduit enables safe and efficient repeated access or removals of tools and tissue



Precise control

- Enables colon shortening, to straighten tortuous anatomy for easier, faster access and improved endoluminal tool usage



Enhanced visualization


- Provides gentle endoluminal anchoring via the balloon which helps reveal hidden tissue planes and polyps

**COMPLEX
EMR**

**ESD & COLONIC
SUTURING**

**COLONIC
EUS**

**CHALLENGING
COLONOSCOPIES**



EXCEPTIONAL GLIDING PERFORMANCE

Lumendi's proprietary hydrophilic coating creates a self-lubricating layer between the sheath and endoscope, making navigating the colon anatomy quicker and easier.

A unique tool for complex endoluminal procedures

DiLumen EZ¹ is a non-sterile accessory designed for immediate availability. Setup requires only activation of the overtube with water or saline, ensuring minimal preparation time.

This device enables rapid and controlled advancement to the right colon, even in anatomically challenging cases.

Once deployed, it offers dependable assistance in enhancing the technical possibilities of endoluminal therapy.



DiLumen EZ¹

Endoluminal Interventional Platform

Product codes:

- | | |
|-----------|--|
| D-1100-5 | DiLumen EZ1 130 cm working length (5 units/box) |
| D-1100-10 | DiLumen EZ1 130 cm working length (10 units/box) |
| D-1101-5 | DiLumen EZ1 103 cm working length (5 units/box) |
| D-1101-10 | DiLumen EZ1 103 cm working length (10 units/box) |

Contact us at
info@lumendi.com to book your personal product demonstration and follow the link to explore clinical results & case reports.

