


DESIGNED TO BENEFIT PATIENTS, PHYSICIANS AND EFFICIENCY



One Hysteroscope for a Range of Solutions

The KARL STORZ 26153EA Hysteroscopy Set offers patients, physicians, and nursing staff a range of advantages in terms of comfort, positive results, and efficiency gains.

This continuous-flow, single-piece hysteroscope sheath is ideal for use with the 2.9 mm telescopes. Exceptional flow dynamics help to ensure that optimum uterine distention is maintained throughout procedures.

This is one hysteroscope that delivers solutions to enhance patient comfort, support positive surgical outcomes, and drive efficiency gains.

STORZ
KARL STORZ – ENDOSKOPE

DESIGNED TO BENEFIT PATIENTS, PHYSICIANS AND EFFICIENCY

One Hysteroscope for a Range of Solutions

Enhance Patient Comfort

- Small 5.4 mm outer diameter makes this sheath ideal for office procedures
- Angled flow parts minimize discomfort to the patient's legs

Supports Positive Physician Outcomes

- 21 cm extended working length provides easy access to tubal ostia
- Sealing Cap prevents fluid leakage and helps maintain a dry environment for physicians
- Compatible with both 12° and 30° telescopes

Drive Efficiency Gains

- Hysteroscope Sheath ensures quick and simple assembly, while helping to minimize the expense of lost or damaged equipment and repairs

Telescopes	
26153 EA	Hysteroscope Sheath , size 5 mm
26153 EAU	Sealing Cap , single use endoscopic seal GYN, sterile, package of 10
26120 BA	HOPKINS® Forward-Oblique Telescope 30° , diameter 2.9 mm, length 30 cm
26020 FA	HOPKINS® Telescope 12° , diameter 2.9 mm, length 30 cm

OPTION 2 - All-In-One System	
TP 100EN	TELE PACK X LED , endoscopic video unit
20 2121 34	TELECAM C-MOUNT One-Chip Camera Head , NTSC
20 2000 43	C-MOUNT Lens , soakable, f = 38 mm
495 NA	Fiber Optic Light Cable , diameter 3.5 mm, length 230 cm

OPTION 1 - Office Video System – Office cart with IV pole, drawer, and T-stand	
20 2330 11	TELECAM DX II Camera Control Unit
20 2121 34	TELECAM C-MOUNT One-Chip Camera Head , NTSC
20 2000 43	C-MOUNT Lens , soakable, f = 38 mm
9702-19	19" Barco LCD Monitor
201315-20	Cold Light Fountain XENON NOVA® 175
9805T-20	Office Cart with IV-Pole , Drawer, and T-Stand
495 NA	Fiber Optic Light Cable , diameter 3.5 mm, length 230 cm

Hand Instruments	
26159 UHW	Biopsy and Grasping Forceps , semirigid, 5 Fr., length 34 cm
26159 BHW	Biopsy Spoon Forceps , semirigid, 5 Fr., length 34 cm
26159 EHW	Scissors , semirigid, blunt, jaws, 5 Fr., length 34 cm
26159 DHW	Punch , semirigid, through-cutting, 5 Fr., length 34 cm
26159 H	HESSELING Tenaculum Grasping Forceps , semirigid, 5 Fr., length 34 cm

Cleaning Accessory	
27650 C	Cleaning Brush , outer diameter 2.5 mm

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com