

SIZE 3.5 PARSONS LARYNGOSCOPE THE RIGHT SIZE FOR YOUR PATIENT



When Size 3 is Too Small & Size 4 is Too Large

Bridge the gap with our size 3.5 PARSONS Laryngoscope

- Provides excellent exposure
- Small enough to fit a wide range of patients
- Large enough to introduce up to a size 6 rigid bronchoscope

STORZ
KARL STORZ—ENDOSKOPE

SIZE 3.5 PARSONS LARYNGOSCOPE THE RIGHT SIZE FOR YOUR PATIENT

When Size 3 is Too Small & Size 4 is Too Large

At times, a size 3 PARSONS Laryngoscope is too small for the best exposure and working space for your instruments, and a size 4 is too large for your patient. Similarly, when a PARSONS Laryngoscope is used to aid the introduction of a size 4.5 rigid bronchoscope, a size 3 PARSONS is too small for the bronchoscope tube and a size 4 is once again too large for your patient. But, our size 3.5 PARSONS Laryngoscope is designed to fit that patient.

“The perfect laryngoscope for large infants to young adolescents.”

— David Parsons, M.D.

Ordering Information	
8576CC	PARSONS operating laryngoscope, size 3.5, length 13 cm



For more information about the size 3.5 PARSONS Laryngoscope, or other KARL STORZ products, contact your local representative or call 800-421-0837

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com